

N95 Guide

For additional masks, contact Processed Stores at 319-356-1784

N95 Model Fit Tested To	PS#	NIOSH Equivalent	PS#
1860	155877	8210	156518
1860S	155878	8110S	157163
1870+	155876	9205+, 9210+	158407, 160124
1804	155879		
1804S	155880		

NOTE: If using NIOSH equivalent and engaging in a procedure that might subject you to high velocity liquid (e.g., an arterial bleed) you must:

- Wear a face shield to protect the NIOSH mask, or
- Wear a medical-grade mask over the NIOSH N95